

BaseSpace Onsite Platform Release Notes

For BaseSpace Onsite v1.2

August 25, 2014

Revision History

Revision	Date	Description of Change
A	8/25/14	Initial Release.

Introduction

These Release Notes detail key changes to software components for BaseSpace between the specific versions listed in the table below.

Software Application	Prior Version	New Version
BaseSpace Onsite	1.1	1.2

BaseSpace Onsite v1.2

NEW FEATURES:

- System and Active Analysis Logging features added to Log tab of BaseSpace Admin screen
 - The admin user can now generate additional logging information to assist with Illumina support team diagnosis

DEFECT REPAIRS:

- Automatic log collection for Illumina support debugging has been fixed.
- Timeout when setting Gateway IP in network settings has been fixed.
- Now suppressing erroneous predictive power failure messaging.